
Laboratory Equipment

in customized premium quality

MADE IN GERMANY

LET'S SHAKE IT

Edmund Bühler GmbH

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Edmund Bühler GmbH

Founded by Edmund Bühler in the university city of Tübingen in 1880, Edmund Bühler GmbH is still owned by the family today.


More than 140 years later, Edmund Bühler GmbH has developed to become a company which operates throughout the world. The production and sales of conventional laboratory devices is traditionally a main focus of the company. A second focal point is the development and production of research equipment, in particular for the quickly growing field of materials science.

In the laboratories of the world, the „Bühler Shaker“ is considered to be a synonym for the efficiency of the company. Competence, experience and the commitment to excellence in quality, which is steadily anchored in the company philosophy, form the basis for a unique program of technically mature products.

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Description

Table-top shaker for high loads and continuous operation

Basic equipment

Basic device with shaking plate and rubber mat, without rack system

Technical Data

Motion SM 30 A:	orbital
Motion SM 30 B:	reciprocating (horizontal)
Motion SM 30 C:	can be switched
	from orbital to reciprocating
Shaking platform:	560 x 400 mm
Max. load:	30 kg
Shaking speed:	15 - 300 rpm, steplessly adjustable
Stroke SM 30 A:	26 mm
Stroke SM 30 B:	30 mm
	46 mm (Option: No. 0052 212)
	50 mm (Option: No. 0052 216)
Stroke SM 30 C:	26 mm
Timer:	0 - 120 min / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 46 W
Heat emission:	approx. 7 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	33 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	42 pieces
Erlenmeyer flasks 250 / 300 ml	20 pieces
Separating funnels 250 ml	6 pieces
Separating funnels 1000 ml	2 pieces
Test tube racks	4 pieces

Further loading possibilities see page 47/48

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**Rack system
Combifix SM A**
No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels.

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2-storey top frame SM
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 323 mm. Distance between the 2 storeys: 282 mm (height sufficient for 1000 ml Erlenmeyer flasks)

Alternatively, the rack systems Combifix SM or the universal tray SM can be mounted on the 2-storey top frame SM.

For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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**Rack system
Combifix SM B**
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v. For fastening horizontal vessels between the clamping strips, or as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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Sliding Plates

Sliding plates with Combifix SM (rack systems Combifix SM, but with additional base plate).

Sliding Plate with Combifix SM A
No. 0051 484

with Combifix SM B
No. 0051 485

with Combifix SM C
No. 0051 486

Mounting set
No. 0052 074
for universal trays, or sliding plates on the 2-storey top frame SM.

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**Rack system
Combifix SM C**
No. 0050 156

Rack system with 2 clamping strips h, 2 cramp strips and 1 spring strip. For fastening separating funnels. Standard version of cramp strip: 5 clamps NS 19 (D20) and 2 clamps NS 29 (D32). With additional strips also suitable for other vessels.

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Universal tray SM
No. 0051 472

Coated tray with drillings for individual fastening of spring clamps or test tube racks. Shaking platform 560 x 400 mm.

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Description

Table-top shaker with programmable shaking tasks for the mixing of difficult media

Basic equipment

Basic device with shaking plate and rubber mat, without rack system. Programmable control with touch display (Technical details of control see page 44)

Technical Data

Motion SM 30 A control:	orbital
Motion SM 30 B control:	reciprocating (horizontal)
Motion SM 30 C control:	can be switched from orbital to reciprocating
Shaking platform:	560 x 400 mm
Max. load:	30 kg
Shaking speed:	15 - 300 rpm, in steps of 5
Stroke SM 30 A control:	26 mm
Stroke SM 30 B control:	30 mm 46 mm (Option: No. 0052 212) 50 mm (Option: No. 0052 216)
Stroke SM 30 C control:	26 mm
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 46 W
Heat emission:	approx. 8 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	33 kg

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USB 2.0 interface

No. 0052 931

for the communication between PC and shaker

RS 232 interface

No. 0052 933

for the communication between robotic system and shaker

Control interface

No. 0052 932

for the integration of the shaker in robotic systems (SPS)

Safety Package*

No. 0060 024

Monitoring of motor and shaker speed

Defined stop*

No. 0060 020

Stop always in the same position (SM 30 A control and SM 30 B control)

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*Installation only possible in our works

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2-storey top frame SM
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 323 mm. Distance between the 2 storeys: 282 mm. (height sufficient for 1000 ml Erlenmeyer flasks)
Alternatively, the rack systems Combifix SM or the universal tray SM can be mounted on the 2-storey top frame SM.
For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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Description

Table-top shaker with heating plate and programmable shaking tasks for mixing of difficult media

Basic equipment

Basic device with heating plate, without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion:	orbital
Shaking platform:	560 x 400 mm
Max. load:	25 kg
Shaking speed:	15 - 300 rpm, in steps of 5
Stroke:	26 mm
Runtime:	programmable / continuous
Temperature range:	up to 120°C
Control accuracy:	approx. +/- 1°C
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Enclosure protection:	IP 21
Power consumption:	max. 450 W
Heat emission:	400 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	40 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	42 pieces
Erlenmeyer flasks 250 / 300 ml	20 pieces
Separating funnels 250 ml	6 pieces
Separating funnels 1000 ml	2 pieces

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**Rack system
Combifix SM B**
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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The below specified shaking speeds against load are approximate values. Depending on specific properties of substances to be shaken these values can differ marginally.

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Description

Programmable shaking tasks for the mixing of difficult media

Basic equipment

Basic device with shaking plate and rubber mat, without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion:	reciprocating (horizontal)
Shaking platform:	560 x 400 mm
Max. load:	30 kg
Shaking speed:	15 - 240 rpm, in steps of 5
Stroke:	50 mm for Gushing Test (No. 0052 216)
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 46 W
Heat emission:	approx. 8 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	33 kg

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2 x Spannleiste h Nr. 0050 400
1 x Spannleiste q Nr. 0050 403

1-storey rack system
No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels. (see picture)

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2-storey top frame
No. 0052 065

Top frame with 2 coated trays with drillings for mounting rack systems Combifix SM or for spring clamps or test tube racks.

Total height: 323 mm;
Distance between the 2 storeys:
282 mm.

Dimensions (w x d x h):
680 x 420 x 480 mm
Weight: 7.7 kg

Alternatively, the universal tray SM can be mounted on the 2-storey top frame SM. For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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3-storey top frame
No. 0039 206

Top frame with 3 coated trays for mounting the rack systems Combifix SM. Total height: 400 mm; Distance between the storeys: 177 mm.

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Maximum Shaking Speed Beer Bottle Shaker SM 30 B

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Nr. / No. 6162 000

Description

Excellent system thanks to reproducible ambient temperatures and useful accessories
Can be combined with all universal shakers or used as independent unit

Basic equipment

Coated metal housing with 3 windows, the front door opens upwards, socket inside, cable feedthrough Ø 12,6 mm at the back

Technical Data

Temperature range:	5°C above ambient to +50°C
Temperature accuracy:	+/- 1% of the set value
Inner dimensions (w x d x h):	660 x 540 x 430 mm (height sufficient for 2000 ml Erlenmeyer flasks)
Lighting:	standard light 15 W / UV light 15 W
Temperature sensor:	PT 100
Air circulation:	240 m ³ /h
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Power consumption:	max. 620 W
External Dimensions (w x d x h):	680 x 610 x 560 mm
Weight:	33 kg

Safety cut-off of UV light and stop of the shaker when the front door is opened if combined with a Bühler Universal Shaker
Standard temperature controller 5°C above ambient to +50°C

Options

Extended temperature range up to +60°C

No. 0052 105

Installation of an additional heater in order to reach temperatures up to max. +60°C. Installation only possible in our works.

Cooling device TH 30

No. 0051 581

Installation of a cooling coil and connectors for an external cooling system (refrigerated circulator or flow-through cooler). Depending on the cooling capacity of the external cooler temperatures of approx. 5°C below room temperature can be reached (300 - 500 W required). Installation only possible in our works.

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Shelf TH 30
No. 0052 123

Additional tray for samples which have to be temperature controlled but need no shaking.
Dimensions: 365 x 558 mm.

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2-storey top frame SM / TH
No. 0052 117

2-storey top frame for all SM shakers. Dimensions adapted to TH 30. For vessels with a max. height of 14 cm (Erlenmeyer flasks 200 / 250 ml).
Dimensions (w x d x h): 600 x 402 x 192 mm.
For fastening the vessels, either universal trays SM or sliding plates with Combifix SM can be used.

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- Intuitive operation via capacitive touch display for easy operation also with gloves
- Non-reflecting LCD colour display 110 x 67 mm with resolution 800 x 480 dot
- Bright LED background lighting (400cd/m²)
- Display of date and real time
- Sensor for the detection of an incubator hood on the shaker
- Manual operation or automatic operation
- Pause function for the shaking tasks
- Logfile function for the recording of interruptions of the shaking process
- Possibility to change the rotational direction for optimal mixing (orbital motion)
- Automatic operation:
 - 10 programs with 10 segments each, freely definable
 - Optical and acoustic signal at the end of a program
 - Possibility to program a defined time for starting the shaking action
 - Programmable restart after power failure
 - Alpha-numeric program names

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Options

USB 2.0 interface (socket type B) for the communication between PC and shaker No. 0052 931
Plug-and-play connection with standard terminal programs.
Complete remote control of the shaker.

Interface RS 232 for the communication between robotic system and shaker No. 0052 933
Serial interface for industrial and robotic systems.
Shaker completely remote controlled.


Control interface for the integration of the shaker in robotic systems (PLC) No. 0052 932
- Relais output: shaker is running / stops
- Relais output: signal alarm / trouble-free operation
- Relais output: ready for operation
- Optocoupler input for external start/stop
- Analog input for setpoint of the speed
- Analog output for actual speed value

Defined Stop*
Predefined stop always in the same position. Especially suited for the use of the device in robotic systems.
(TiMix 5 control, KS 15 A control, SM 30 A control, SM 30 B control, VKS 75 A control; KS 15 B control and VKS 75 B control on request)

Safety Package*
The speed of both motor and shaker is continuously monitored. In case of deviations, e.g. due to a defective drive belt, a warning message is displayed and an acoustic signal is audible.

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Clamping strip h

Holder with sponge rubber cover for fastening upright vessels.

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Clamping strip v

Holder with sponge rubber cover in combination with clamping strip h.

For fastening horizontal or high vessels.

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Cramp strip

For fastening horizontal separating funnels.

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Spring strip

For securing the stoppers of separating funnels or other horizontal vessels.

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For universal trays. The sizes are related to Erlenmeyer flasks, but are also suitable for round-bottom flasks, beakers, etc.

10 ml	25 ml	50 ml	100 ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml
0009 642	0009 643	0009 644	0009 645	0009 646	0009 647	0009 648	0009 649	0009 653	0009 652

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To be fastened either with clamping strips on rack systems or with a hinged foot on universal trays. Dimensions 265 x 100 mm each.

With the hinged foot it is possible to vary the angle of inclination of the test tubes.

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Abmessung 289 x 100 mm.

Required for universal trays. Size fitting for all test tube racks. Dimensions 289 x 100 mm.

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¹⁾ with 5 clamping strips h (0052 115) additional ²⁾ with 4 clamping strips h (0052 115) additional ³⁾ with 3 clamping strips h (0052 115) additional
⁴⁾ with 2 clamping strips h (0052 115) additional ⁵⁾ with 2 clamping strips h (0050 118) additional ⁶⁾ with 1 clamping strip h (0050 118) additional
⁷⁾ with 3 clamping strips h (0050 118) additional ⁸⁾ with 3 clamping strips h (0050 400) additional ⁹⁾ with 1 clamping strip h (0050 400) additional
¹⁰⁾ with 5 clamping strips h (0050 387) additional ¹¹⁾ with 1 clamping strip h (0050 387) additional

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¹²⁾ with 1 clamp strip (0050 206) and 1 spring strip (0050 207) additional ¹³⁾ Number of pieces depends on shape and dimensions of the separating funnels

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Examples of customized solutions

With our comprehensive laboratory device program we can offer coherent answers for most demands which arise from practical laboratory work. We believe that our job is even more. A special challenge are the very particular areas of use for which the standard is not enough. Here we design and produce customized solutions. Customer-oriented, quality conscious and tailor made. Some examples of our customized products:

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Rotationsschüttler RS 10
Rotation Shaker RS 10

If undelivered please return to below address

Book Post

Scientific Research Instruments Company Private Limited

42-43, 2nd & 3rd Floor, 1st Cross, Gubbalala, Bengaluru - 560 061

Telephone : 080 4757 2577 | 080 4302 5791

Sales

info@srico-labworld.com

Mobile : +91 9900674407

Service

service@srico-labworld.com

Mobile : +91 9900055879

Bhubaneswar | Hyderabad | New Delhi | Navi Mumbai | Vadodara

Ahmedabad | Chandigarh | Chennai | Goa | Guwahati | Kolkata | Lucknow | Pune | Thiruvananthapuram | Visakhapatnam

www.srico-labworld.com